

CASE ID

ARCH

☐ Max ☐ Mand ☐ Both

SURGERY TYPE

☐ Guided ☐ Freehand ☐ Hybrid

PROVISIONAL PLAN

☐ Immediate ☐ Early ☐ Delayed

RECORDS CHECKLIST

- ☐ Approved setup reference (diagnostic wax-up or new setup)
- ☐ DICOM dataset complete and readable
- ☐ STL files submitted and labeled consistently
- ☐ Implant system documented
- ☐ Scan-body brand + exact reference number(s) documented
- ☐ Soft-tissue capture submitted
- ☐ Opposing record submitted
- ☐ Bite relationship submitted
- ☐ Vertical dimension reference documented
- ☐ Smile line and midline references documented
- ☐ Occlusal scheme goal noted
- ☐ Provisional reference (if present) submitted
- ☐ Post-op deviations communicated to Burbank Dental Lab

A missing scan-body reference number can freeze a full-arch case even when every scan file is already submitted — the geometry that transfers implant position is reference-specific, and one brand can have several designs that aren't interchangeable.