

BIOHORIZONS Implant	Torque (Ncm)	Torque Settings	Driver & Head
Internal 3.5, 4.5, 5.7	30 Ncm	0.050" (1.25 mm) Hex	
Tapered Internal 3.5, 4.5, 5.7	30 Ncm	0.050" (1.25 mm) Hex	
External Implant 3.5	30 Ncm	0.050" (1.25 mm) Hex	
External Implant 4.0, 5.0	30 Ncm	0.050" (1.25 mm) Hex	

BIOMET 3i Implant	Torque (Ncm)	Torque Settings	Driver & Head
Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	20 Ncm	0.048" Hex	
Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	20 Ncm	0.048" Hex	
Certain™ XP 4/5, 5/6	20 Ncm	0.048" Hex	
Miniplant™ 3.25	20 Ncm	0.048" Hex	
Osseotite™ 3.75, 4.0, 5.0, 6.0	20 Ncm	0.048" Hex	
Osseotite™ XP 4/5, 5/6	20 Ncm	0.048" Hex	

DENTSPLY Implant	Torque (Ncm)	Torque Settings	Driver & Head
ANKYLOS® C/X - C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)	15 Ncm	1 mm Hex	
ANKYLOS® C/X - /X 3.5, 4.5, 5.5, 7.0 (indexed)	15 Ncm	1 mm Hex	
OsseoSpeed™ EV 3.0 - Green	25 Ncm	0.050" Hex	
OsseoSpeed™ EV 3.6 - Purple	25 Ncm	0.050" Hex	
OsseoSpeed™ EV 4.2 - Yellow	25 Ncm	0.050" Hex	
OsseoSpeed™ EV 4.8 - Blue	25 Ncm	0.050" Hex	
OsseoSpeed™ EV 5.4 - Brown	25 Ncm	0.050" Hex	
OsseoSpeed™ Profile EV 4.2 - Yellow	25 Ncm	0.050" Hex	
OsseoSpeed™ Profile EV 4.8 - Blue	25 Ncm	0.050" Hex	
OsseoSpeed™ / OsseoSpeed™ TX EV 3.0 - Yellow	25 Ncm	0.050" Hex	
OsseoSpeed™ / OsseoSpeed™ TX 3.5, 4.0 - Aqua	15 Ncm	0.050" Hex	
OsseoSpeed™ / OsseoSpeed™ TX 4.5, 5.0 - Lilac	20 Ncm	0.050" Hex	
OsseoSpeed™ 4.0 - Blue (Previous)	25 Ncm	0.050" Hex	

DENTSPLY (Cont.) Implant	Torque (Ncm)	Torque Settings	Driver & Head
OsseoSpeed™ TX Profile 4.5, 5.0	25 Ncm	0.050" Hex	
XIVE® S 3.0	24 Ncm	1.22 mm Hex	
XIVE® S 3.4, 3.8, 4.5, 5.5	24 Ncm	1.22 mm Hex	
FRIALIT™ 3.4, 3.8, 4.5, 5.5	24 Ncm	1.22 mm Hex	

NOBEL BIOCARE® Implant	Torque (Ncm)	Torque Settings	Driver & Head
Angulated Screw Channel	35 Ncm	Nobel Biocare® Omnigrip	
NobelActive™ 3.0	15 Ncm	Nobel Biocare® Unigrip	
NobelActive™/NobelReplace™ Conical Connection NP 3.5	35 Ncm	Nobel Biocare® Unigrip	
NobelActive™/NobelReplace™ Conical Connection RP 4.3, 5.0	35 Ncm	Nobel Biocare® Unigrip	
NobelReplace™ NP 3.5	35 Ncm	Nobel Biocare® Unigrip	
NobelReplace™ RP 4.3	35 Ncm	Nobel Biocare® Unigrip	
NobelReplace™ WP 5.0	35 Ncm	Nobel Biocare® Unigrip	
NobelReplace™ 6.0	35 Ncm	Nobel Biocare® Unigrip	
Brånemark System™ NP 3.3	35 Ncm	Nobel Biocare® Unigrip	
Brånemark System™ RP 3.75 4.0	35 Ncm	Nobel Biocare® Unigrip	
Brånemark System™ WP 5.0, 5.5	35 Ncm	Nobel Biocare® Unigrip	

STRAUMANN® Implant	Torque (Ncm)	Torque Settings	Driver & Head
Bone Level 3.3 NC (Narrow CrossFit™)	35 Ncm	SCS Screwdriver	
Bone Level 4.1, 4.8 RC (Regular CrossFit™)	35 Ncm	SCS Screwdriver	
Standard/Standard Plus 4.8 RN (Regular Neck)	35 Ncm	SCS Screwdriver	
Standard/Standard Plus 4.8 WN (Wide Neck)	35 Ncm	SCS Screwdriver	
Standard Plus 3.3 NN (Narrow Neck)	35 Ncm	SCS Screwdriver	

ZIMMER Implant	Torque (Ncm)	Torque Settings	Driver & Head
Tapered Screw-Vent 3.5, 4.5, 5.7 / 3.5, 4.5, 5.7 friction fit	30 Ncm	1.25 mm Hex	
Screw-Vent 3.3, 3.7, 4.5 / Screw-Vent 3.3, 3.7, 4.5 friction fit	35 Ncm	1.25 mm Hex	
SwissPlus™ 4.8	35 Ncm	1.25 mm Hex	